STABILITY TESTS FOR 99MTC-DTPA RADIOPHARMACEUTICAL

PROF. ASOC. BRUNILDA DACI PROF. DR. ELIDA BYLYKU PROF. ASOC. KOZETA TUSHE PROF. ASOC. DRITAN PRIFTI

RAP.24

INSTITUTE OF APPLIED NUCLEAR PHYSICS UNIVERSITY OF TIRANA

INTERNATIONAL CONFERENCE ON RADIATION APPLICATIONS

- Applications for marketing authorization in respect of radiopharmaceuticals should be accompanied, as in case of all medicinal products, by the particulars and documents referred to in Directives 65/65EEC and 75/31/EEC as amended and in the Annex of Directive 75/318 EEC as amended.
- Most radiopharmaceuticals are used for the purpose of medical diagnosis. They contain only small amounts of active substances with a radionuclide attached to them to allow scintigraphic imaging or measurement of biodistribution.

- Radiopharmaceuticals have changing composition with time, associated with radioactive decay. The physical half-life of the radionuclide is often so short (like ^{99m}Tc). In these cases, the final preparation must be done immediately before administration to the patients.
- For all radiopharmaceuticals, the shelf life of the product as supplied by the manufacturer should be specified and justified, as should a shelf life after reconstitution where applicable, considering radiochemical and radionuclide degradation products.

- The aim of this study was determination of stability of ^{99m}Tc-DTPA (Diethylenetriaminepentaacetic acid) radiopharmaceuticals used for brain and kidney scanning.
- Stability is defined by analyzing radiochemical purity of the ^{99m}Tc-DTPA radiopharmaceuticals in different times after reconstitution (labelling of cold DTPA-kit with ^{99m}Tceluate) and production (shelf-life of cold DTPA-kit).


- Chromatographic methods of assaying radiopharmaceuticals are used for determining the radiochemical purity of ^{99m}Tc-DTPA radiopharmaceuticals.
- The system used was ITLC-SG/Acetone and ITLC-SG/Saline to determine the sum of radiochemical impurities (free and hydrolyzed ^{99m}Tc). In conclusion of many experiments, the results show that ^{99m}Tc DTPA radiopharmaceutical was stable up to 2 hours after its reconstitution and 1 year after cold kit production.

DTPA - RADIOPHARMACEUTICAL

SUMARY DATA

STRUCTURAL FORMULA


- Chemical Name:
 Diethylene Triamine Penta Aceticacid (DTPA)
- <u>- Technetium tc-99m DTPA</u> is a Radioactive Diagnostic Agent.
- The mechanism of action of <u>technetium tc-</u> <u>99m DTPA</u> is as a Radiopharmaceutical Activity.
- A technetium diagnostic aid used in renal function determination



PROCEDURE

- Eluation of ⁹⁹Mo-^{99m}Tc Generator
- Labelling of cold kit
- Developing of chromatography in different time after labelling and production
- Scanning of the radiochromatograms
- Analyzing the results

ELUATION OF GENERATOR




LABELLING OF COLD KIT

Content (%) of the compounds in cold kit

Compounds	Version A	Version B
DTPA	<mark>5.00 mg</mark>	<mark>9 mg</mark>
SnCl ₂ x2H ₂ O	<mark>0.25 mg</mark>	<mark>0.6 mg</mark>
Ascorbic acid	<mark>0.15 mg</mark>	<mark>0.5 mg</mark>
NaCl	1	<mark>20 mg</mark>







STABILITY OF LABELLED KIT

Stability in-vitro means for how long the labelled kit is stable, which means that radiochemical purity more than 90%.

 Determination of radiochemical purity 15 min, 3 hours and 24 hours after labelling.

The results ate presented in the table below

RADIOCHEMICAL PURITY

- Radiochemical purity is the fraction of total activity in the desired radiochemical form present in the radiopharmaceutical.
- In the radiopharmacy, it is the presence of the undesired radiochemical impurities that is checked.
- These impurities are due to decomposition of the radiopharmaceutical caused by solvent, temperature, light of radiolysis or labeling of a chemical impurity with the same radionuclide.

Radiochemical purity (% 99m Tc-R/F) = 100 – (% 99m TcO₄⁻ - % 99m TcO₂)

CHROMATOGRAPHIC SYSTEMS

INSTANT THIN LAYER CHROMATOGRAPHY METHOD

Stationary Phase:


ITLC-SG

Mobile Phase:

Acetone and NaCl

DEVELOPMENT OF CHROMATOGRAPHY











RADIOCHEMICAL PURITY

Activity curve (PC) Acetone 85%




Activity curve (PC) 0.9%NaCl %



RADIOCHEMICAL PURITY

Activity (%) of cut parts

DTPA (Acetone)		DTPA (0.9%NaCl)		
Origin	97.00 (%)	Origin	2.61 (%)	
Middle	1.92 (%)	Middle	1.81 (%)	
Front	1.18 (%)	Front	95.58 (%)	



STABILITY OF LABELLED KIT

In - Vitro Stability of kit prepared/labelled (two versions)

Time after	991	m Tc – DTPA (Ver	sion A)	^{99m} Tc – DTPA (Version B)			
labelling ^{99m} Tc free	^{99m} Tc coloid	^{99m} Tc complex	^{99m} Tc free	^{99m} Tc coloid	^{99m} Tc complex		
15 min	0.12	0.35	99.53	0.22	0.43	99.35	
1 hour	0.41	0.34	99.25	0.29	0.40	99.31	
2 hours	0.50	0.40	99.10	0.43	0.41	99.16	

SHELF-LIFE

Shelf-life of DTPA (3 months after production)

Time after labelling	^{99m} Tc – DTPA (Version A)			^{99m} Tc – DTPA (Version B)		
	^{99m} Tc coloid	^{99m} Tc free	^{99m} Tc complex	^{99m} Tc coloid	99mTc free	^{99m} Tc complex
30 min	3.10	2.80	94.10	1.79	0.09	98.12
2 hours	2.95	3.10	93.95	2.04	1.07	96.89

Shelf-life of DTPA (1 year after production)

Time after labelling	^{99m} Tc – DTPA (Version A)			^{99m} Tc – DTPA (Version B)		
	^{99m} Tc coloid	^{99m} Tc free	^{99m} Tc complex	^{99m} Tc coloid	^{99m} Tc free	^{99m} Tc complex
15 min	3.5	2.64	93.86	1.93	0.16	97.91
2 hours	1.90	4.12	93.98	1.30	0.27	98.43

CONCLUSIONS

- . 99m TC DTPA Radiopharmaceutical produced in sterile and freeze dried form fulfill all the requirements for stability
- -` The results show that ^{99m}Tc DTPA radiopharmaceutical was stable
 - up to 2 hours after its reconstitution and
 - 1 year after cold kit production in freeze dried form
 - ITLC method used to determine radiochemical purity is accurate, simple and fast

THANK YOU